

STUDY OF SURFACE PASSIVATIONS ON HPGe DETECTORS

Serena Riccetto
University of Perugia – INFN Perugia

2nd Position Sensitive Germanium Detectors (PSeGe) Workshop – Milano (Italy), September 2017



STUDY OF SURFACE PASSIVATIONS ON HPGe DETECTORS


Outlines

- Characterization by γ -spectroscopy of a planar HPGe detector prototype
 - Study of the passivated surface with low energy γ -ray
 - Comparison between different passivation processes:
 - Methanol passivation
 - Sulfide termination
 - Low H: hydride termination in 10% concentration acid
 - High H: hydride termination in 50% concentration acid
 - Hyper H: same as High H but with an additional etching
-
- Previous work
- This work
- NEW

New experimental set-up

Several HPGe detector passivations have been realized in the LNL's Materials Laboratory and a new experimental setup has been constructed to allow a better study of their properties.

- Always the same planar HPGe n-type crystal (21 mm height - 39 mm diameter)
- Reduced lateral thickness of the mounting cap
- Lead collimator on a micromotional device




Measurements


After every passivation procedure, the detector was mounted in the special cryostat, cooled and subjected to a complete characterization with γ sources.

Two different kinds of measurements:

- determination of detector bulk properties by putting different standard calibration γ sources in front of the end cap;
- scanning of the lateral intrinsic surface by irradiating from the side of the crystal with an ^{241}Am collimated source on a micrometrical motion device that moves with 0.6 mm steps




Bulk properties



Determination of:

- the best response working voltage or depletion voltage V_d , defined as the reverse bias voltage required to extend the depletion region through the full thickness of the semiconductor bulk
- the energy resolution for the 1.33 MeV ^{60}Co peak as a function of voltage
- the detector efficiency as a function of the radiation energy


Bulk properties




Depletion voltage V_d

- The depletion voltage is not the same for the different passivations: 2000V for the sulfide, 1700V for the methanol and 1500V for the two hydrides
- The counting rate at the plateau is different for the sulfide: it was proved to be related to a smaller active volume

- ^{60}Co source, in front of the detector
- Integrals of the peak at 1.33 MeV, normalized to the acquisition live times




Bulk properties




Depletion voltage V_d

- The depletion voltage is essentially the same for these three passivation methods, as well as the final plateau value

- ^{60}Co source, in front of the detector
- Integrals of the peak at 1.33 MeV, normalized to the acquisition live times




Bulk properties




Energy resolution

- The best resolution does not depend on the passivation, it is essentially the same for all the methods [1.74-1.79 keV]

- ^{60}Co source, in front of the detector
- Resolution (FWHM) of the peak at 1.33 MeV



Bulk properties




Energy resolution


➤ The best resolution is essentially the same for all these three methods, as well as the voltage at which those resolution are reached

Passiv.	Depletion Voltage [V]	Best resolution Voltage [V]	Best resolution [keV]
Low H	1500	1500	1.74
High H	1500	1500	1.75
Hyper H	1500	1500	1.77

- ^{60}Co source, in front of the detector
- Resolution (FWHM) of the peak at 1.33 MeV




Bulk properties



Efficiency

- The efficiency shows not to depend on the H-type passivation process and it is comparable with the expected values and with the results of the other passivations


- Different standard sources, in front of the detector




Passivation properties

Lateral scanning with collimated source for the:

- determination of the dead layer below the lateral surface induced by the passivation, estimated starting from the counting rate and applying the Beer-Lambert law for γ rays absorption
- study of the Am spectra as a function of the distance between the collimated γ beam and the detector contacts, for each passivation treatment at its best resolution voltage and in some cases at different voltages




Passivation properties




- S passivation is the most used in microelectronics because it's stable, but it has the larger dead layer



- ^{241}Am collimated source, from the side of the detector
- Integrals of the full-energy peak at 59.54 keV, normalized to the acquisition live times




Passivation properties




The counting rate in:

- Low-H passivation decreases moving from the n^+ to p^+ contact, giving rise to a slightly n-type surface;
- High-H passivation increases from n^+ to p^+ contact, giving rise to a slightly p-type surface;
- Hyper-H passivation remains constant along all the scan, giving rise to an ideal passivation.

- ^{241}Am collimated source, from the side of the detector
- Integrals of the full-energy peak at 59.54 keV, normalized to the acquisition live times




Passivation properties




The counting rate in:

- Low-H passivation decreases moving from the n^+ to p^+ contact, giving rise to a slightly n-type surface;
- High-H passivation increases from n^+ to p^+ contact, giving rise to a slightly p-type surface;
- Hyper-H passivation remains constant along all the scan, giving rise to an ideal passivation.

- Dead layer thickness determination from the counting rates applying the Beer-Lambert law for gamma absorption




Passivation properties




Background and resolution


- An increase of the background contribution to the spectra with the dead layer thickness is observed for all the different passivations



- ^{241}Am collimated source, from the side of the detector




Passivation properties




Background and resolution


- An increase of the background contribution to the spectra with the dead layer thickness is observed for all the different passivations



- ^{241}Am collimated source, from the side of the detector




Passivation properties






Background and resolution

- An increase of the background contribution to the spectra with the dead layer thickness is observed for all the different passivations



- ^{241}Am collimated source, from the side of the detector

HYPER H



Conclusions

- The efficiency, the depletion voltage and the energy resolution of the detector are comparable for the three H-terminated passivations
- The estimated dead layers induced below the passivated surfaces are thinner and more homogeneous than the ones of common commercial detectors, but they have presented different trends: the Hyper H termination is the most uniform

Collaborators

INFN-LNL

INFN-LNL and University of Padova

INFN-LNL and University of Verona

INFN-LNL and University of Trento

INFN-PG and University of Camerino

INFN-PG and University of Perugia

IKP Cologne

CNR-IMM Bologna

D.R. Napoli, W. Raniero (ENSAR2)

D. De Salvador, G. Maggioni, S. Carturan,
E. Napolitani, V. Boldrini, F. Sgarbossa

G. Mariotto

G. Della Mea

N. Pinto

S. Ricetto

J. Eberth

R. Nipoti, F. Mancarella, M. Bellettato



THANK YOU FOR THE ATTENTION