


Laser Spectroscopy in an MR-ToF Device

Anthony Roitman


McGill University (Canada), CERN (France/Switzerland)

L

- Major area of nuclear research: stability at shell closures
 - ▶ Reflected in many observables, such as binding energy or charge radius




kink in charge radii at shell closures




X. Yang et al., Prog. Part. Nuc. Phys. 129, 104005 (2023)

kink in charge radii at shell closures





disappearance of shell closure at N=20: island of inversion



X. Yang et al., Prog. Part. Nuc. Phys. 129, 104005 (2023)

- S. J. Novario et al., Phys. Rev. C 102, 051303 (2020)
- D. T. Yordanov et al., PRL, 108:042504 (2012)






X. Yang et al., Prog. Part. Nuc. Phys. 129, 104005 (2023)

- S. J. Novario et al., Phys. Rev. C 102, 051303 (2020) D. T. Yordanov et al., PRL, 108:042504 (2012)
- N = 20 shell closure disappears for magnesium: charge radii for 33,34 Mg would provide a powerful benchmark for nuclear theory


ISOLDE facility




ISOLDE facility



Challenge: yields down to 100 ions / second and half lives down to 40 ms



Challenge: yields down to 100 ions / second and half lives down to 40 ms



Need highly sensitive laser spectroscopy techniques to probe charge radii of exotic nuclei


Laser Spectroscopy in Nuclear Physics

By probing an atom's electronic hyperfine structure, we can determine the properties of its nucleus, such as:

nuclear spin





electromagnetic moments



charge radii









• Collinear geometry minimizes Doppler broadening:


$$\delta \nu \propto \frac{\delta E}{\sqrt{F}}$$



ullet Limitation: Measurement time $ullet\sim 5~\mu s$, but $ullet T_{1/2}>10~ms$




Solution?



Solution?

Multi-Ion Reflection Apparatus for CLS (MIRACLS)


Multi-Reflection Time-of-Flight (MR-ToF) device increases effective beampath to "recycle" ions



• signal-to-noise ratio improvement: $\frac{S}{N} = \frac{S_0}{N_0} \sqrt{r}$


Multi-Ion Reflection Apparatus for CLS (MIRACLS)

Multi-Reflection Time-of-Flight (MR-ToF) device increases effective beampath to "recycle" ions




- signal-to-noise ratio improvement: $\frac{S}{N} = \frac{S_0}{N_0} \sqrt{r}$
- More exotic radionuclides with low production yields can be probed


• Single-passage mode (experimental data):




• Multi-reflection improvement (experimental data):




• Multi-reflection improvement (experimental data):




• Multi-reflection improvement (experimental data):




Beam energy: large source of uncertainty in CLS



Beam energy: large source of uncertainty in CLS




Beam energy: large source of uncertainty in CLS



$$\Rightarrow \nu_0 = \sqrt{\nu_c \cdot \nu_{\text{a}}},$$
 independent of beam energy!

Beam energy: large source of uncertainty in CLS



Incoming laser
(lasers switched



 $\Rightarrow \nu_0 = \sqrt{\nu_c \cdot \nu_a}$, independent of beam energy!

Nuclear Charge Radii of ^{7,9,10}Be and the One-Neutron Halo Nucleus ¹¹Be


W. Nörtershäuser^{1,2}, D. Tiedemann², M. Žáková², Z. Andjelkovic², K. Blaum³, M. L. Bissell⁴, R. Cazan², G. W. F. Drake⁵, and Ch. Geppert^{1,6} et al.


MIRACLS method

▶ A short animation


Status January 2022







Highest energy MR-ToF device ever built! (> 10 kV)




Highest energy MR-ToF device ever built! (>10~kV) Advantages for mass separation: see F. M. Maier et al., NIMA A., 1048, (2023)

November 2024: Mg Beamtime


Measured our first physics case: ³⁴Mg





November 2024: Mg Beamtime


Measured our first physics case: ³⁴Mg



34Mg measurement with down to 10 ions in MR-ToF device per bunch

November 2024: Mg Beamtime


Measured our first physics case: ³⁴Mg




- 34Mg measurement with down to 10 ions in MR-ToF device per bunch
- ³³Mg measurement also performed with repumping scheme.

Sensitivity Limit


Performed measurements with down to 0.7 ions / measurement cycle in MR-ToF device




Goal: Extend previous measurements by the COLLAPS experiment to 98,99 Cd (N = 50 shell closure)



Results:




Results:



• Also 10 ions in MR-ToF per bunch, but with much lower laser power.

Results:



- Also 10 ions in MR-ToF per bunch, but with much lower laser power.
- Pump & probe scheme used for ^{99,101,103}Cd

Summary

 MIRACLS has demonstrated its effectiveness as a new technique for CLS

Summary

- MIRACLS has demonstrated its effectiveness as a new technique for CLS
- Magnesium: charge radii measured for even-even isotopes ^{24–34}Mg and odd-even ³³Mg

Summary

- MIRACLS has demonstrated its effectiveness as a new technique for CLS
- Magnesium: charge radii measured for even-even isotopes ^{24–34}Mg and odd-even ³³Mg
- Cadmium: charge radii measured for even-even isotopes ^{98–116}Cd, charge radii and moments for odd-even ^{99,101,103}Cd

Acknowledgements



Collaboration:











Funding:










MIRACLS team and participants:

L. Croquette, H. Heylen, E. Leistenschneider, S. Lechner, F. Maier, L. Nies, P. Plattner, M. Rosenbusch, F. Wienholtz, M. Vilen, R. Wolf, F. Buchinger, W. Nörtershäuser, L. Schweikhard, S. Malbrunot-Ettenauer O. Ahmad, T. Fabritz, P. Giesel, R. Hernandez, H. Heylen, J. Hughes, F. Koehler, K. Koenig, D. Lange, L. Lalanne, T. Lellinger, E. Matthews, A. Mcglone, K. Mohr, J. Palmes, V. Repo, L. Rodriguez, C. Shweiger, J. Spahn, J. Warbinek, J. Wilson, Z. Yue, C. Zambrano

Odd-even scheme for magnesium

=


• Repurposing Lasers: One for repumping, one for scanning



Odd-even scheme for magnesium

=

• Repurposing Lasers: One for repumping, one for scanning



Odd-even scheme for cadmium

