Search for the QCD Axion with ADMX

The Axion Quest- Recontres du Vietnam August 5, 2024

Nick Du on behalf of ADMX Postdoctoral Scholar





Outline

- Axion Dark Matter
- Axion Haloscopes
- ADMX Run 1D
- Plans for Future Searches





Axion Dark Matter

- The nature of dark matter is one of the largest mysteries in physics
- Peccei-Quinn solution to the Strong CP problem predicts the axion
 - Axions have properties that make them compelling as a candidate for dark matter
- Axion dark matter would solve two huge mysteries in physics at once
 - This is very compelling!



TODAY



Photo: https://sanpedrosupermarket.com/product/axion-lime-dish-soap/



Axion Landscape



Adapted from Ciaran O'Hare



Axion Dark Matter Searches: Haloscopes



Adapted from Ciaran O'Hare





Axion Haloscopes





Scanning with Axion Haloscopes

- Axion mass is not known *a priori*
 - Axion haloscopes must be tunable across a wide frequency range
- Tuning rods set the resonant frequency of the cavity and are used to sweep the cavity resonant frequency target frequency range
 - Power from cavity is sampled then the cavity resonance is tuned and the procedure is repeated







Scan Rate for Cavity Axion Haloscopes

$$\frac{df}{dt} \approx 323 \frac{\text{MHz}}{\text{year}} \left\{ \left(\frac{g_{\gamma}}{0.36}\right)^4 \left(\frac{f}{1 \text{ GHz}}\right)^2 \left(\frac{\rho_0}{0.45 \frac{GeV}{cc}}\right)^2 \right\} \cdot \left\{ \left(\frac{3.5}{SNR}\right)^2 \left(\frac{B_0}{7.6 \text{ }T}\right)^4 \left(\frac{V}{136l}\right)^2 \left(\frac{Q_L}{30,000}\right) \left(\frac{C_{lmn}}{0.4}\right)^2 \left(\frac{0.35 \text{ }K}{T_{sys}}\right)^2 \right\}$$

Parameters fixed by Nature

- g_{γ} Dimensionless Coupling constant
- f Axion frequency
- ρ_0 Dark matter halo energy density

Experimental Parameters

- SNR Signal-to-noise ratio
- B_0 External magnetic field
- *V* Cavity volume
- Q_L Cavity quality factor
- C_{lmn} Cavity form factor
- T_{sys} System noise temperature



Axion Dark Matter Experiment (ADMX)





ADMX Run 1D

- Run 1D
 - Began January 2024
 - $-\,$ Search for axion in the 1-1.4 GHz frequency range
- Resonator design
 - Same tube but with larger tuning rod to reach a higher frequency
- Receiver design
 - Tunable current-pumped JPA enables quantum-limited noise across entire range
 - Anticipating transition to a flux-pumped JPA by Washington University







Data Taking Run with ADMX

- Data taking operations are controlled by an automated script
- Multiple measurements are required at each tuning step to monitor status and optimize the noise temperature of the experiment







Synthetic Axion Injections

- Synthetically generated axion waveforms are injected into the cavity via a weakly coupled antenna
- Enables testing of RF receiver and candidate identification methodology
 - Blinded axion injections
- We regularly inject synthetic axions into our system to ensure a robust test of our

identification methodology





Step 3:



Vetoing Candidate Signals

LLNL-PRES-867546

- Candidates are always expected due to blind injections, noise fluctuations, RFI
 - Candidate veto procedure ensures reject non-axion signals





Upcoming Limits on the Axion with ADMX Run 1C Extend



LAWRENCE Livermore National Laboratory



Probing for Higher Mass Axions

- The target mass of your axion search sets the length scale of your resonant cavity
- As resonant frequency of the cavity goes up
 - Volume decreases as $V \sim 1/f^3$
 - Quality factor decreases as $Q \sim 1/f^{2/3}$
 - Noise power increases at $T_{amp} {\sim} f$
- To maintain an adequate scan rate need new developments (Multiple cavities, Stronger magnets, Higher Q cavities, etc.)

$$\frac{df}{dt} \approx 323 \frac{\text{MHz}}{\text{year}} \left\{ \left(\frac{g_{\gamma}}{0.36}\right)^4 \left(\frac{f}{1 \text{ GHz}}\right) \left(\frac{\rho_0}{0.45 \frac{\text{GeV}}{\text{cc}}}\right)^2 \right\} \cdot \left\{ \left(\frac{3.5}{\text{SNR}}\right)^2 \left(\frac{B_0}{7.6 \text{ T}}\right)^4 \left(\frac{V}{136l}\right)^2 \left(\frac{Q_L}{30,000}\right) \left(\frac{C_{lmn}}{0.4}\right)^2 \left(\frac{0.35 \text{ K}}{\text{T}_{sys}}\right)^2 \right\}$$

$$\text{Lower loss} \text{Below SQL noise (?)}$$

$$\text{magnets} \text{Larger cavities to} \text{maintain volume} \text{improve Q factor}$$



ADMX Sidecar: High Frequency Testbed

- Higher frequency cavity mounted above the main cavity
 - Testbed for new resonator designs and technology
- Sidecar is also an axion search!





Bartram et. al., Rev. Sci. Instrum. 94. 044703 (2023)



ADMX Sidecar: Hybrid Superconducting Cavities

- Hybrid cavity will utilize superconducting tuning rod with copper cavity
 - Superconducting tuning rod will reduce signal losses from tuning rod
 - Expecting a 30% increase in quality factor
- Using Nb₃Sn thin-film super conductor
 - Demonstrated to maintain superconductivity in high magnetic fields
 - Provided by S. Posen at FNAL







ADMX 2A: Multi-Cavity Searches

- Multiple cavities combined in a phasematched method scan scale the SNR of an axion by the number of cavities
- Run 2A: 1.4-2 GHz
 - Utilize 4 co-tuned cavities
 - Tuning will be done with piezoelectric motors for faster frequency matching
- Extended Frequency Range: 2-4 GHz
 - 18 cell cavity array
 - To be run in 9.4 T MRI magnet at Fermilab





ADMX-EFR: Cryostat Design







ADMX-EFR: Magnet Installation

- Magnet has been installed at Fermilab
- Currently undergoing diagnostics before it will undergo ramping







ADMX-EFR Cavity Prototype

- ADMX EFR will consist of an array of 18-cavities
 - Cavities are 1 m long, 128 mm diameter
 - Cavities will be horizontally mounted inside the magnet
 - Tuning rod armature acts as a counterweight for the tuning rod
- Different diameter tuning rods will change the frequency range
 - 32 mm rods: 2.0-3.1 GHz
 - 54 mm rods: 3.1-3.9 GHz









ADMX-EFR Tuning Impedance Filter

- Stepped impedance filter on the axles of tuning rods prevent coupling of power out of cavity by tuning rod axles
 - Enables fabrication of metal tuning rod axles, instead of classical dielectric
 - Easier fabrication
 - Improved thermalization
 - Reduced dielectric losses





ADMX-EFR: Tuning System

- Each tuning rod is rotated by a rotary piezoelectric motor coupled to the cavity via a gear
 - Angular resolution on the rotary motors is $1 \ \mu deg$
- A linear piezoelectric motor will adjust the insertion depth of a dipole antenna to control coupling to the cavity
- We can vary coupling and rod position with fine control





ADMX-EFR: Two-Cavity System Cooldown

- Cryogenic testing of the EFR prototype is underway at LLNL!
- Cavities were cooled to 100 mK
- Two-cavity testbed has been useful to demonstrate modelocking and tuning in a cryogenic environment







Summary

- Axion dark matter is exciting field!
- ADMX is currently searching for axion dark matter in the 1.0-1.4 GHz frequency regime
- Multicavity arrays will improve on the sensitivity for higher frequency axion searches





Acknowledgements

This work was supported by the U.S. Department of Energy through Grants No DE-SC0009800, No. DE-SC0009723, No. DE-SC0010296, No. DE-SC0010280, No. DE-SC0011665, No. DEFG02-97ER41029, No. DE-FG02-96ER40956, No. DEAC52-07NA27344, No. DE-C03-76SF00098 and No. DE-SC0017987. Fermilab is a U.S. Department of Energy, Office of Science, HEP User Facility. Fermilab is managed by Fermi Research Alliance, LLC (FRA), acting under Contract No. DE-AC02-07CH11359. Additional support was provided by the Heising-Simons Foundation and by the Lawrence Livermore National Laboratory and Pacific Northwest National Laboratory LDRD offices.



Questions





Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor Lawrence Livermore National Security, LLC, nor any of their employees makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or Lawrence Livermore National Security, LLC. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government purposes.

