

# An overview of optical benches for mirror metrology and Coating Thermal Noise (CTN) measurement



Raffaele Flaminio, François Frappez, Enzo Galizzi, Romain Gouaty,  
Chinmay Malik, Alain Masserot, Sylvain Petit, Balazs Sashalmi,  
**Antoine Syx**, Edwige Tournefier

ET France Workshop - 01/04/26

# Birefringence characterization of Sapphire substrates

Sapphire is a potential candidate for ET but is a **birefringent** material.


- Development of a bench to measure birefringence.
  - ◆ Adapted to ET sized mirrors.

Principle of operation:

- A **10.5 mW laser** source (1064 nm or 1550 nm) is **polarized** using a Glan-laser calcite polarizer.
- A Photoelastic Modulator (PEM) **modulates the polarization** state of the beam at 50 kHz.
- The **beam** interacts with the **sample** where it experiences a **phase delay** due to **birefringence**.
- A second polarizer converts the **phase delay** to an **intensity modulation**.
- A lock-in **amplifier** detects this modulation to **isolate the signal** at the PEM modulation frequency.
- Measure **retardation magnitude** ( $\Delta n \times d$ ) and **fast-axis orientation** ( $\rho$ ).



Figure 1: Birefringence bench setup ongoing



Credit: Severin Nadji for LMA (for ET Core Optics meeting in March 2026)



# Study micro-scale absorbing defects of mirror coatings

Detect and characterize point absorbers in mirror coatings.

→ Shared bench with the birefringence measurement (see previous slide).

Principle of operation:

- A **200 W laser** (1064 nm) is directed towards the **mirror**.
    - Point absorbers create thermal gradients.
  - A **1 mW probe laser** (532 nm) **scans** the **mirror surface**.
  - A Shack-Hartmann **sensor** captures **distorted wavefronts**.
  - Comparison with thermal camera to **identify point absorbers**
- First measurements allowed to assess uniform absorptions above 100 ppm and detect point absorbers.



Next measurement: test on two Virgo's IM spare mirrors.

# Absorption measurement for crystalline coatings


ANR MICRONG



This bench was updated from the previous state of measuring absorption on amorphous samples to be able to measure 200 mm crystalline wafers.

Principle of operation:


- A **high power laser** (1064 nm) is directed towards the **mirror**.
  - Substrate absorption induces thermal gradients/refractive index gradient.
- A **probe laser** (1300 nm) is pointed at an angle on the **same mirror area**.
  - Wavelength is **transmitted** by the coating.
  - Beam is **deflected** due to the thermal gradient.
- The probe beam's deflection is detected.



Cartographie Ø150 mm - pas de 5 mm AC1940

Measure AlGaAs 200 mm:

7.4 +/- 7.9 ppm




# CTN measurement using Cyclic Interferometry

Time-differential measurement of CTN without locking a cavity.

Principle of operation:

- A **laser** (NIR) is injected into a **Sagnac interferometer**.
  - A 3x3 coupler splits it **clockwise (CW)** and **counterclockwise (CCW)**.
- The beam travels through the **delay loop** (L=4 km).
  - A **time delay** is induced.
- Each beam interacts with the mirror.
- The beam **recombines** at the coupler.
  - **Interference** signal sensitive to **phase differences**.
- **CTN** induces **phase difference** between beams.
  - Detected by a balanced photodiode.



# Development of an optical bench for the measurement of CTN using a resonant cavity

ANR MICRONG

# Measuring the CTN (MIT) - The cavity

Direct CTN measurement using a laser beam in a Fabry-Pérot cavity.

How to isolate the mirror sample's CTN from the cavity mirrors' CTN?

- Use of a 3-mirrors folded cavity with the sample mirror in the center.
- The effect of the CTN on the beam increases as the beam size on the mirror decreases.
  - A large beam size on the input and output mirrors is required.
  - And a smaller beam size on the mirror sample.




Figure from [Phys. Rev. D 95, 022001](#)

# Measuring the CTN (MIT) - Multimode readout

How to isolate the CTN from parasitic noises?

S. Gras and al.  
[Phys. Rev. D 95, 022001](#)  
[Phys. Rev. D 98, 122001](#)

- Many parasitic noises will dominate over the CTN.
  - The laser frequency and power noise.
  - The cavity length variations (due to vibrations).
- Use of co-resonating transverse modes in the cavity : TEM00 (fundamental), TEM02 et TEM20 (high order modes “**HOM**”).
  - Dominant parasitic noises common to each modes.
  - **HOM** TEM02 and TEM20 scan different areas of the sample.
  - Extraction of the CTN = measuring the frequency difference between the modes 02 and 20.




Figure from [Phys. Rev. D 95, 022001](#)





# Study of a test cavity - OMC

We test the setup with a test cavity.

An Output Mode Cleaner (OMC):

- Monolithic silica block.
  - Less sensitive to vibrations than the MIT cavity.
- 4 reflective surfaces = mirrors.
  - We measure the total thermal noise of the cavity (including thermo-refractive noise which is expected to be dominant in this case).




# Generation of the beams


## Pre-fiber part

Laser Nd-YaG (1064 nm) with 1 W of power split in 3 beams:

- Beam 1 (TEM00) phase modulated with fibered Electro-Optic Modulator (EOM) at 18.75 MHz.
- Beam 2a (TEM02) frequency shifted using Acousto-Optic modulator (AOM) at 200 MHz.
- Beam 2b (TEM20) frequency shifted using AOM at 208 MHz.
- Each beam is injected in a fiber.



Current state




Layout






# Intensity stabilization of beams 2a & 2b

- **Beams 2a & 2b** are output from the fibers.
- For each beam, a fraction of the power is detected by a **photodiode**.
- The signal is used to stabilize the intensity of the beam by acting on the **AOMs amplitudes**.
- Intensity Stabilization Servo (**ISS**) implemented analogically (UGF @ 40 kHz).




# Frequency lock of beams 2a & 2b to the cavity's TEM02/20

- Using a **quadrant photodiode (QPD)**, 2a & 2b beams are frequency locked to the **TEM02/TEM20** of the cavity by acting on the **AOMs frequencies**.
- Control loop implemented digitally (UGF @ 1.3 kHz).
  - To be implemented analogically to improve the bandwidth.




# Measuring the CTN using the Beat Note signal of beams 2a & 2b


- Beams 2a & 2b interfere on **beat note(s)** photodiodes.
- The frequency difference between the two beams produce a beat note signal at **~ 8MHz**.
- The **CTN** is extracted by demodulating this signal with an external source at **~8 MHz** and measuring frequency fluctuations.




**Lock at 70% of the modes peak = best sensitivity to measure frequency fluctuations**



# Current state of the setup at LAPP



# Goals of the MICRONG experiment at LAPP

- **In-air** measurement of the OMC.
  - Tested with low finesse OMC (280\*) but signal not sensitive enough to measure CTN.
  - **Now** testing with higher finesse OMC (10000\*).
- Building a **3 mirror folded cavity** (work in progress). 
- Moving under **vacuum** of the **cavity** for more test.
- Moving under **vacuum** the **detection + cavity** part (post-fibers).

