

Master Project « Targeted Radiation Therapies »


(2026-2029)

- **CNRS Nuclei & Particles
IN2P3**

**Scientific coordinators : Rachel DELORME (LPSC, Grenoble),
Anne-Marie FRELIN-LABALME (GANIL, Caen), Patrick VERNET
(LPCA, Clermont-Ferrand)**


Tech. Coord: Amelia Maia LEITE (IJCLab, Orsay)

rachel.delorme@lpsc.in2p3.fr




Context: Targeted Radiation Therapies


➤ **Therapeutic strategies for aggressive or systemic cancers**
= molecular targeting of cancer cells




Targeted Radionuclide Therapy (TRT)



Accelerator-based Boron neutron capture therapy (BNCT)




NP-enhanced RT (radiosensitizers)



- **Internal irradiation** (full or secondary reaction)
- Need for a “**vector**” with specific targeting
- **Short range** of particles, potential **high-LET**
- **Dosimetry challenges:** imaging quantification, high heterogeneity, control → *adapted num & instru tools*
- **Neutron production / detection**
- **Biological mechanisms** may differ from EBRT

	α-TRT	BNCT	NP
Radionuclide/particle	^{223}Ra , ^{225}Ac , $^{212/213}\text{Bi}$, ^{211}At , ^{212}Pb ...	$^{10}\text{B}/^{11}\text{B}^*$	e- (PE, Auger)
Energies α (& ^7Li) or e-	5-9 MeV	0.8-1.7 MeV	0-100 keV
Range α (& ^7Li) or e-	40 – 100 μm	5 – 9 μm (<cell)	0-100 μm
LET (keV/μm)	60 – 100	≥ 200	0.5-20

MP Targeted Radiation Therapies: organization





NUCLÉAIRE
& PARTICULES

MP Targeted Radiation Therapies: organization

& Partners

(academic, clinical, industrial)



IN2P3 contributors

~8 FTE
researchers

LPSC

S. Chabod
D. Dauvergne
R. Delorme*
M.L. Gallin-Martel
D. Santos

Subatech

N. Beaupère
T. Bossis (Postdoc)
D. Thers

GANIL

A.-M. Frelin Labalme*

IP2I

M. Alcocer Avila
M. Beuve
L. Cousin (PostDoc)
C. Della Rosa (PhD)
C. Rodriguez-Lafrasse
C. Strubbia (PhD)
E. Testa

LPCA

I. Balandier
G. Farge
L. Maigne
G. Montarou
P. Vernet*

IPHC

F. Bertile

IJCLab

L. Hardouin (PhD)
P. Lanièce
L. Ménard
Q. Mouchard
O. Seksek
M.A. Verdier

LPC-Caen

M. Rousseau

~2 FTE
Tech & Ing

LPSC

J.F. Muraz
M. Reynaud
C. HOARAU
J. F. Muraz
N. Sauzet
O. Guillaudin
J. Marpaud

IP2I

G. Alfonse
Céline Malesys
G. Garde
Y. Zoccarato
R. Della Negra
BE Mécanique

LPCA

P.-E. Vert
M. Lefebvre
G. Rivrais
D. Grimbichler

IJCLab

A. Maia Leite*
D. Crepin
A. Said
A. Semsoum
A. Bouafia
R. Leplat

*coordinators

WP1: Deciphering biological responses from molecules to cells in targeted radiotherapies

Characterize and quantify the biological specificity of targeted RT efficacy, through detailed radiobiological studies using radionuclides, NP or high-LET particles (relevant for α -IRT and BNCT).


T1.1 Functional impact of α -emitters and nanoparticles on molecular targets

Molecular and functional analysis

- Online analysis of molecular and functional damage following DNA metabolizing enzymes exposure with α -emitters
- Mitochondrial / metabolism impact of dose enhancement due to metal NPs.

Follow-up of irradiation impact on cells (*instrumental dev. in WP3*)

- In Clermont-Ferrand: Mini-X prototype of an epifluorescence microscope coupled with a mini-irradiator (X-rays)
- In GANIL: scintillator-based imaging setup + microscope for α -IRT exp.



Hibernation: bio-inspired approach to radiosensitize human cells

- Deciphering specific mechanisms involved in radiosensitization
- Characterization of bear serum specific compounds



WP1: Deciphering biological responses from molecules to cells in targeted radiotherapies

Characterize and quantify the biological specificity of targeted RT efficacy, through detailed radiobiological studies using radionuclides, NP or high-LET particles (relevant for α -TRT and BNCT).

T1.1 Functional impact of α -emitters and nanoparticles on molecular targets

T1.2: Biological data production to feed and benchmark models

Presentation Alicia
Garnier (mardi 31 Juin)

In vitro biological effect of β -emitting isotopes (*IRHydrobrain*)


Strategy for intraoperative TRT in glioblastoma with radiolabeled chitosan hydrogel with ^{90}Y or ^{177}Lu .

- Radiobiological effects of hydrogel will be studied on F98-spheroid model in terms of DNA damage (*gH2AX foci*), cell depth, cell migration, ROS measurements.

In vitro biological effect of α -emitting isotopes (*AlphaBioDose*, *OPTADOM*)

Application to breast cancer brain metastasis treated with α -IRT (and other cell lines)

- Clonogenic cell survival and sub-cellular damage (DNA, mitochondria, membrane) measurements, on various cell lines, with α radionuclides (^{223}Ra , ^{212}Pb ...) thanks to dosimetrically instrumented experiment
- and with low-energy ^4He and ^7Li ions on SiLab and BioALTO platforms



WP2: Multiscale modeling of physical dose and biological response

2.1: Modeling at the subcellular, cellular, and multicellular scales:

→ Develop **advanced modeling tools** to reproduce TRT irradiation dose distributions with heterogeneity complexity and **predict biological damage** at **cellular**, and **multicellular** scales.

Main obj. 2026:
3 Master internship
Poster at PTCOG

Past and current developments:

- **NanOx Biophysical model:** predict cell survival probability from physical dose distribution with mix of **nano- and micro-dosimetry** and **consideration of oxydative stress**.

→ Recently **adapted** from hadrontherapy to **low-energy ion TRTs***

- **Application to α -TRT:**



Calculation pipeline: **CPOP - Geant4 - NanOx**

to study impact of **intra-cell** and **inter-cell heterogeneities** of sources on **cell survival** and **tumor control probability** in spheroids.

→ **effect++ of inter-cell activity fluctuation****

- **Preliminar application to BNCT:**

→ **Huge impact of intracellular heterogeneities**




*Alcocer-avila et al., 2024, Med. Phys. 51(12):9358.

**Levrage et al., 2025, Med. Phys. 52:e17917

Main perspectives (~4y):

1. **Benchmark of calculated cell survival** to biological in vitro experiments (in TRT & BNCT).

2. **Impact study of realistic cell morphologies:** microscopy-reconstructed numerical phantoms



3. **Intra-cell & intra-tumoral heterogeneity impact study for BNCT** (and GdNCT).

4. **NanOx extension to new cancerous cell lines** (U87, SQ20, MDA...) and **extra-nuclear damage**.

5. **GATE bio-dose actor developments** applied to α IRT and BNCT treatments

WP2 – Including LQD Cross Sections in G4-DNA


Presentation Loic
COUSIN (mardi 31 Juin)

LQD Code : Monte-Carlo simulation of Ions transport in liquid water



In collaboration with CIMAP (Caen), Rosario university (Argentina) and IP2I (Lyon)

Objective : *Including LQD cross sections (CS) in a new Physic List of G4-DNA.*

LQD : Electron Cross Sections




LQD : Ion Cross Sections



For Ions (H, He, C, ...):

- Ionisation :
 - CDW-EIS Model :
 - Simple Diff and Double Diff CS : (E_i vs E_{sec}, Angle)
- Multiple Ionisations
- Outlooks : Uses the new G4 LQD Physic List within the collaboration



Ref : Gervais, M. Beuve, G.H. Olivera, M.E. Galassi,
Numerical simulation of multiple ionization and high LET effects in liquid water radiolysis,
<https://doi.org/10.1016/j.radphyschem.2005.09.015>.


WP2: Multiscale modeling of physical dose and biological response

2.2: Preclinical dose and biological effect modeling

➔ Aims at developing **modern dosimetry workflow** based on **multi-modal biodistribution imaging**, **Monte Carlo simulation** and **statistical approaches**, to optimize experiments and interpret preclinical results.

Development of GATE modeling tools to support innovative β -IRT strategy for intraoperative GBM (IRHydrobrain):

GATE 10 used to support **in vitro** & **in vivo** experiments:




3D reconstruction (spheroid, mouse, rabbit and human) from microscopic or CT imaging.
Geometries in GATE.

1. Calculate **dose to targets** with **rat digital twin** (*in vivo correlation*)
2. Implementation of a **spheroid digital twin** (*in vitro exp.*)
 - 2.1 implementation of a **Chemistry actor**
 - 2.2: **ROS quantification** and comparisons
 - 2.3: **DNA damage** and comparisons

Application to breast cancer early brain metastasis treatment with α -IRT (AlphaBioDose, OPTADOM):


1. **Activity biodistribution determination** of antibodies in **mice** using **multimodal information**: PET imaging and autoradiography, vs anatomical MRI and histology.



2. **Dose rate and dose distribution calculation** using GATE
3. **3D biological dose maps calculations** (*from T2.1 TCP models*) and *in vivo* response comparison.
4. **Learning by statistical methods** for dose and bio-dose reconstruction

WP1/2/3 – 3D dosimetry in Targeted Radionuclide Therapy (TRT)


- Dosimetry workflow applied to a clinical case (S. Lechien PhD thesis)*



*Application
to preclinical
studies*

Experiment
performed in
2026 with ^{223}Ra



- Assessment and integration of biological effect of alpha particles (A. Doudard PhD)*





Scanning electron microscopy micrograph of micro-metastatic clusters and illustration of α and β^- particles

A emitters:

- Short range
- High LET
- Hypoxia sensitivity



Main achievement: developed a **Si detector** for precise dosimetry in *in vitro* α -emitter experiments, using a deconvolution of energy spectra to reconstruct radionuclide z-position distribution in the wheel.



Major article: Doudard A et al., *Application of a new spectral deconvolution method for in vitro dosimetry in assessment of targeted alpha therapy.* Med Phys. 2023;50:3762–3772. <https://doi.org/10.1002/mp.16279>

WP3: Optimization and control of dose delivery

3.1: Optimization and characterization of neutron fields in BNCT

→ Aims at proposing an **optimal set of accelerator-based BNCT systems**, including the study of **innovative targets** for the production of intense epithermal neutron field and their test system, the **design of optimized moderators**. and **instrumental developments for neutron field detection** and microdosimetry.

Main obj. 2026:

Experiment + ICNCT conf

in Oct 2026 Buenos Aires

Past and current developments:

1. Innovative targets for intense epithermal neutron:

- Developed a rotating target on a graphite wheel + ^9Be (or ^{13}C) deposit
- e- thermal test beamline (3 kW/cm²)




2. Neutron field detection

- NFM gas detector with ^{10}B -foil for **field monitoring**
- MIMAC-FastN **neutron spectrometer in epithermal and fast ranges** (*angular distributions meas. legnaro*)



3. Optimal moderator (BSA) design:

- Developed **new topology optimization (TopOpt)** method → deepest treatment depth achievable in the AB-BNCT community (*~10 cm vs std ref ~7.6cm*)




Main perspectives

1. Dev. of a rotating ^{13}C (d (1.5 MeV,n)) target for the Buenos-Aires AB-NCT accel.
2. BSA design: **Optimal Moderator design** for the ^{13}C (d (1.5 MeV,n)) reaction + **for metrological epithermal neutron** Cadarache (LMDN-ASNR)
3. **Neutron spectro**: charac. The TANDAR field
4. Test MIMAC as **microdosimeter** with Tissue equivalent gaz
5. **NFM detector validation** on metrological ASNR neutron high flux & at GANIL
6. Ultimate goal: propose **optimal AB-BNCT facility design**

WP3: Optimization and control of dose delivery

3.2: Dose delivery control in *in vitro* and preclinical radiobiology

➔ Aims at developing dedicated irradiation facilities and detectors enabling accurate and quantitative assessment of biological experiments, in full complementarity of WP1.




BioALTO irradiation platform (~4y perspectives):


Aims to provide **multiple low- and medium-energy ion beams**, with irradiation quality **compatible with radiobiology studies**.

1. **Upgrade RG device** to improve ALTO interoperability
2. **Develop new beam diagnostics and dosimetry tools:**


Diamond counter



Si-based
μdosemeter



Air
fluorescence
Profiler



3. **Creation and equipping a cell culture room @ALTO**
4. **Dosimetric and radiobiological validation**
5. **Feasibility studies for UHDR irradiations**

Presentation
Louis Hardouin

Mini-X for online functional analysis: Coupling between a microscopy platform and a mini-irradiator (X-rays)


- enables **real-time visualization** of the **early cellular responses to irradiation**.
- **Online dosimetric characterization system** under development




In vitro dosimetry of α-emitters:

Si-based dosimeter designed for use in cell culture incubator ➔ determine real dose received to cells during exp.

- dose information at the scale of the culture well
- **new system to enable microdosimetry** measurements at the scale of **individual cells**. Combines **scintillator-based imaging** (*α tracks visu*) setup with a **microscope** (*damage follow-up*)
- Complementary to mini-X



Frelin-Labalme et al., Med Phys. 2020;47(3):1317.



Preclinical dosimetry in IRT:

Development of a **customized biomimetic mouse phantom** instrumented with **miniature scintillating fibre dosimeters**

➔ for IRT dosimetry evaluation.





NUCLÉAIRE
& PARTICULES




Hope first
dosimetric exp. By
end 2026 / begin
2027

Lauréat SESAME 2025 Région Île de France:

- Jouvence ligne BioALTO
- Jouvence ALTO
- Financement moniteur PEPITES dédié
- Financement microdosimetre dédié

Plan CAO de la nouvelle ligne
BioALTO en cours de montage:



Équipes:



Contact: Amélia Maia Leite, Quentin Mouchard
bialto@ijclab.in2p3.fr

WP3: Optimization and control of dose delivery




3.3: Clinical dose delivery control in internal radiotherapy (personalized dosimetry)

➔ Aims at developing advanced imaging technologies to enable the individual determination of absorbed doses to organs-at-risk and target regions.

THIDOS :


Clinical prototype of a mobile γ camera of 10x10cm field of view, to control the dose delivered during ^{131}I (β 364 keV) treatment for thyroid diseases.

2025: First clinical evaluation at IUCT.
And in 2 other clinical centers



Perspectives (~1-3y):

1. Adaptation for dose-based treatment planning of thyroid diseases using ^{123}I (160 keV, Cochin)
2. Open use to medium-energy γ -emitters (200-400 keV, ^{177}Lu or ^{225}Ac)
3. Develop a DL-enhanced Compton camera to allow dosimetric monitoring of TAT (Eg $\geq$ 400 keV). ^{149}Tb , ^{213}Bi , ^{225}Ac or ^{212}Pb



XEMIS2:

First liquid Xenon medical imaging system for small mammals

Installed @ Nantes Université Hospital


Dedicated to multimodal imaging (high-energy SPECT, PET, 3X)

Allow image acquisition with reduced activity

Main Perspective for TRT:

Explore it's use for IRT monitoring through images of gamma rays.

➔ feasibility study with ^{212}Pb



AIDER: European project (HORIZONEURATOM, Spain, Italy...):

Development of a Compton Camera for TRT (α -IRT and BNCT):


@Lyon (CREATIS-IP2i) : Camera optimization using MC simulations + image reconstruction

WP3 – Clinical dose delivery control in TRT

- *First clinical evaluation of the mobile THIDOS gamma camera:*

Finalization of the **first clinical evaluation of the THIDOS mobile camera** conducted on **20 patients with benign thyroid diseases** (Claudius Régaud Institute, Toulouse).

Preliminary results indicate **unprecedented image quality** at energies of approximately 400 keV.



Presentation
Marc-Antoine
Verdier (Lundi
29 Juin)

Principales réalisations 2025/2026

- Fin 1er essai en ^{131}I sur patients en mars 26 avec ICR, Toulouse
→ analyse dosi en cours
- Design d'un nouveau collimateur pour le ^{177}Lu , essai prévu avec CHUV Lausanne

Recent articles related to the THIDOS camera:

Bossis et al. (2024). A high resolution portable gamma-camera for estimation of absorbed dose in molecular radiotherapy. IEEE TRPMS, 8(5):556–570


Bossis, T. et al. (2023). Optimized reconstruction of the position of interaction in high-performances γ -cameras. NIM-A, 1048:167907




NUCLÉAIRE
& PARTICULES

Some detailed projects & backup

WP2 – biophysical modeling at microscale in TAT and BNCT



α -emitters
or ^{10}B




Calculation pipeline:

Spheroid modeling = **CPOP**.

Particle transportation = **Geant4**.

Cell survival & TCP conversion from deposited energies = **NanOx**.

TAT: Tumor control probability (TCP) as a function of cell internalization and inter-cell fluctuations of α emitters



• Main Results (V. Levraque PhD thesis & PICTURE project):

Calculation pipeline: Objective to quantify the uncertainties on therapeutic effect predictions due to intracellular and intratumoral heterogeneities

TAT → Lognormal α distribution at tumor scale have a major impact on efficacy. Uniform hypothesis may cause overestimation of TAT efficacy

2 major publications in Medical Physics :


- **Levraque et al**, Impact of intracellular radionuclide distribution in a Monte Carlo biophysical 3D multi-cellular model for Targeted Alpha Therapy. *Medical Physics* **2025**; 52:e17917, <https://doi.org/10.1002/mp.17917>
- **Alcocer et al.**, Biophysical modeling of low-energy ion irradiations with NanOx. *Med Phys.* **2024**;1-14. <https://doi.org/10.1002/mp.17407>

BNCT → Larger impact of cell internalization of ^{10}B on TCP, compared to TAT.

Nuc → Mem = + 500% on TCP

(preliminar results only for Boron dose component)

BNCT: TCP as a function of cell internalization of ^{10}B



→ Perspectives to benchmark model with *in vitro* biological data for TAT & BNCT.

Design of AB-BNCT Beam Shaping Assemblies at LPSC


AB-NCT Challenges:

- To bring $10^9 \text{ n}_{\text{ep}}/(\text{cm}^2.\text{sec})$ to the tumor (30' irradiation time) using a **compact and safe facility** installed at a hospital
- To build a target coping with 30-50 kW on 10 cm^2 (3-5 kW/cm²)
- Shaping and characterizing the neutron field getting out the targets and the moderators = **Beam Shapping Assemblies (BSA)**


Beam Shaping Assemblies (BSAs) are key components of the therapy, positioned between the source and the patient, **to shape and optimize the neutron field used for treatment.** Their design addresses a **complex constrained optimization problem**, until now always solved by parametric optimization.

At LPSC, we have introduced a **new design method, called topology optimization (TopOpt)**, developed internally. This method can compute with little or no approximation the best possible structure of a particle propagator. See S. Chabod et al., Phys. Med. Biol. 70(3) (2025) 035008.

Topology optimized BSAs generate advanced neutron fields, $\varphi(\underline{r}, E, \underline{\Omega})$, that combine neutron directionality ($\underline{\Omega}$), space expansion ($\underline{r}$) and energy distribution (E) to control tumors, generating **treatment depths 30% greater than the best value achieved so far by the global AB-BNCT community.**



Neutron Therapeutics®: scheme of the AB-BNCT treatment facility at Helsinki University Hospital. Left: treatment room, **with the BSA and patient bed.** Right: accelerator and beam line. The line ends with the target, inserted into the BSA.





TopOpt introduces new and unique components into BSA design: neutron tunnels, filters and ring collimators that mimic the effects of multi-beam BNCT treatments.

BNCT literature	Source	TD (cm)
Riley et al. 2003	MIT Reactor	7.8
Kononov et al. 2004	⁷ Li(p(2.3 MeV),n)	7.0
Kim et al. 2009	⁷ Li(p(2.5 MeV),n)	6.0
Inoue et al. 2014	⁹ Be(p(8.0 MeV),n)	5.0
Minsky et al. 2014	⁷ Li(p(2.3 MeV),n)	7.0
Capoulat et al. 2014	⁹ Be(d(1.45 MeV),n)	7.8
Capoulat et al. 2017	¹³ C(d(1.45 MeV),n)	7.8
Torres-Sánchez et al. 2021	⁷ Li(p(2.1 MeV),n)	7.85
Bae et al. 2022	⁹ Be(p(10 MeV),n)	7.2
Rong et al. 2024	⁷ Li(p(2.8 MeV),n)	6.5
Verdera et al. 2024	⁴⁵ Sc(p(2.918 MeV),n)	7
Our work	⁷Li(p(2.5MeV),n)	10.3

Compilation of BSA treatment depth (TD) performance data obtained by the global AB-BNCT community over the past 20 years. TD had progressed by only 0.05 cm in 20 years.


MIMAC-FastN = a fast neutron directional detector to characterize the neutron fields spectrally



700 mbar
He/CO₂ (5%)


Angular distribution
for an AB-BNCT nuclear reaction

Spectrum measured at **0 deg**



INFN LNL
(Legnaro - Italy)

Spectrum measured at **60 deg**



M.E. Capoulat, N.Sauzet *et al.*

« Neutron spectrometry of the $^9\text{Be}(d(1.45 \text{ MeV}),n)^{10}\text{B}$ reaction for accelerator-based BNCT » NIM B, vol. 445, pp. 57-62, 2019

NFM (Neutron Flux Monitor) Detector


Giving the neutron flux (up to 10^{12}n/s) as a function of time, with an evaluation of the dead

time




Multi-layer deposition PAPVD (plasma-assisted)

LPSC – Alexandre BES, Ana LACOSTE

Including a natural Boron layer well measured (thickness, surface and homogeneity) !!



Distribution of the first measurement of the 3D-tracks produced by an epithermal and thermal neutron field at AMANDE (Cadarache) in 2025, in The frame of the **CHEMIN** project (ASNR – CNRS)



Context: TRT Dosimetry

- Radionuclide RT is increasing a lot in clinical practice.
But sub-optimal fixed activity is applied for curative treatment:

- **Dosimetry may guide personalization** to achieve durable complete response by adapting activity / cycle
- **Requires biokinetic imaging and dose calculation models**



➔ Need to simplify dosimetry protocols with the improvement of imaging systems and reconstruction algorithms



- Additional difficulties with high-LET/short range emitters:
Difficulty to observe clear “dose-effect relationship”.

- **High heterogeneity** at micro/tissue scale, often unknown → *homogeneous hypothesis may lead to errors*
- **High (and variable) Relative Biological Efficacy (RBE).**
- **Complex cellular or systemic biological response.**


- We can work on some of the challenges:

- Efficient gamma cameras for personalized dosimetry (> 300 keV)
- Multiscale modeling tools to account for heterogeneities in bio-dose prediction
- Dedicated controlled radiobiological experiments (in vitro / in vivo) to deepen understanding of TRT biological effects

See Ali Ouadi talk for Radionuclides




Garin et al., SIRT using personalized dosimetry for locally advanced hepatocellular carcinoma (HCC) patients. *Lancet*, 6(1),2021




Example of heterogeneous activity distribution in tumor 7min, 7h, and 21h after injection in TAT. Bäck T et al. *J Nucl Med*. 2010;51(10):1616-23.

Context: Accelerator-Based NCT



➤ New era for clinical trials since compact accel. neutron sources (CANS) allow performing BNCT in hospital.

- About **26 new AB-BNCT** facilities world-wide. Clinical trial started in **Finland 2025**, clinical routine for recurrent **H&N cancers** in **Japan**.
- Companies proposing integrated systems (*e.g. neutron therapeutics, TAE life science...*)
- Improved patient care and **beam spectra** for **treating deeper tumors**



Ex. of NuBeam® system (ex. installed in Finland, Japan, UK...)

➤ Some challenges:

- **Production targets:** High power (30-75 kW) degrades targets. Residual radioactivity with ^7Li , needs frequent change. → **Design optimal targets.**
- **Online monitoring:** target aging (neutron flux $\searrow$) and delivered dose → **Develop n & γ (480 keV) detectors.**
- **Beam shaping assembly** → **Find optimal moderation to maximize penetration in tissue** with high-enough intensity.
- **Lack of standardization** → **Neutron field spectral and fluence characterization**
- **Complex bio-dosimetry**, several components.
BNCT biologically-weighted dose based on fix RBE(or CBE) factors:

$$D_W = w_f D_f + w_t D_t + w_\gamma D_\gamma + w_B D_B$$

→ **Models to account for dose micro-heterogeneities & RBE variation**

