

The Target Laboratory at INFN – LNS and the participation in the EURO-LABS project

**Vasilis Soukeras^{1,2}, Manuela Cavallaro¹, Antonio Massara¹, Daniela Calvo³, Antonio Comite⁴, Maria Fisichella^{1,5}, Mauro Giovannini⁴, Martina Ursino¹
for the WP2.5.2 Collaboration**

¹Istituto Nazionale di Fisica Nucleare – Laboratori Nazionali del Sud (INFN – LNS), Catania, Italy

²University of Catania, Catania, Italy

³INFN– Sezione di Torino, Torino, Italy

⁴University of Genova and INFN – Sezione di Genova, Genova, Italy


⁵*Present address:* Grand Accelérateur National d'Ions Lourds (GANIL), Caen, France

***Second meeting on “Targets for Nuclear Physics”
GSI Helmholtzzentrum für Schwerionenforschung, Darmstadt, Germany***

Target: a key – ingredient for a Nuclear Physics experiment

□ A target is an essential ingredient for a nuclear physics experiment !



- ✓ Usually, a beam impinge on a target material, leading to a variety of reaction ejectiles to be measured



Data taking
& analysis

Synergy between involved Researchers and Technicians is imperative !

Targets for Nuclear Physics experiments: solid, gas or liquid ?

- Solid targets are the most common case  **Why ?**
 - ❖ Most elements and materials useful as targets are solid in standard conditions of pressure and temperature.
 - ❖ Much simpler setup
- Gas and liquid targets may be used in some cases  **Why ?**
 - ❖ Need of an active target
 - ❖ Specific properties of some targets

□ Here we focus:

- on solid targets developed and characterised at INFN –LNS 
- on the participation in the EURO-LABS project 

in collaboration with GANIL
R. Rahali et al.



The target lab at INFN – LNS: fabrication

- The target laboratory at INFN –LNS in Catania has over thirty years of experience in target preparation for nuclear physics experiments and interdisciplinary physics
- Users are not only local scientists but also researchers from other universities and laboratories worldwide

A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

Equipment for target production

- ✓ L300 Thermal Bell Jar (TBJ) evaporator
- ✓ L560 Leybold-Heraeus evaporator
- ✓ UNIVEX 400 Leybold evaporator
- ✓ Cold rolling mill

The target lab at INFN – LNS: fabrication

- The target laboratory at INFN –LNS in Catania has over thirty years of experience in target preparation for nuclear physics experiments and interdisciplinary physics


A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

Equipment for target production

- ✓ L300 Thermal Bell Jar (TBJ) evaporator
- ✓ L560 Leybold-Heraeus evaporator
- ✓ UNIVEX 400 Leybold evaporator
- ✓ Cold rolling mill

L300 TBJ evaporator

- 2 resistive sources
- A probe to monitor the backing temperature
- A quartz crystal micro balances



The target lab at INFN – LNS: fabrication

- The target laboratory at INFN –LNS in Catania has over thirty years of experience in target preparation for nuclear physics experiments and interdisciplinary physics


A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

Equipment for target production

- ✓ L300 Thermal Bell Jar (TBJ) evaporator
- ✓ L560 Leybold-Heraeus evaporator
- ✓ UNIVEX 400 Leybold evaporator
- ✓ Cold rolling mill

L560 Leybold-Heraeus

- Evaporation by
 - ✓ e-beam heating source
 - ✓ resistive heating source
- Quartz crystal micro balance
- Halogen heating elements to fix temperature



The target lab at INFN – LNS: fabrication

- The target laboratory at INFN –LNS in Catania has over thirty years of experience in target preparation for nuclear physics experiments and interdisciplinary physics

A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

Equipment for target production

- ✓ L300 Thermal Bell Jar (TBJ) evaporator
- ✓ L560 Leybold-Heraeus evaporator
- ✓ UNIVEX 400 Leybold evaporator
- ✓ Cold rolling mill



UNIVEX 400 Leybold evaporator

- Commissioning tests almost completed!

New

The target lab at INFN – LNS: fabrication

- The target laboratory at INFN –LNS in Catania has over thirty years of experience in target preparation for nuclear physics experiments and interdisciplinary physics

A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

Equipment for target production

- ✓ L300 Thermal Bell Jar (TBJ) evaporator
- ✓ L560 Leybold-Heraeus evaporator
- ✓ UNIVEX 400 Leybold evaporator
- ✓ Cold rolling mill



**Cold
rolling
mill**


LNS Target lab website with some examples of targets produced in the lab:
<https://www.lns.infn.it/en/targets.html>

The target lab at INFN – LNS: characterisation


Equipment for target characterisation

A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

- ✓ Chamber for Alpha-particle Characterisation of target Thickness and Uniformity by Scanning (CACTUS)



- ✓ CACTUS allows for the characterisation of the targets in terms of thickness, local and global non – uniformity.
- Technique: Alpha – particle spectroscopy (APS)







- Chamber diameter $\sim 1\text{m}$
- Host different types of target frame
- Scan different regions of the target surface with a **high precision** (1 mm) thanks to a rotational system and a video camera

The target lab at INFN – LNS: characterisation

A. Massara et al., EPJ Web of Conferences 285, 06003 (2023)

Equipment for target characterisation

✓ Energy calibration based on 2 α -sources (^{148}Gd & ^{241}Am) covering a broad energy range



Target characterisation at INFN – LNS

Why characterisation ?

- ✓ measure the average thickness
- ✓ determine the non – uniformity of the samples

Characterisation is essential towards the selection of the most suitable manufacturing procedure !


How ?

- Alpha Particle Spectroscopy (APS) technique using the CACTUS facility
- Comparison with simulations to determine the **local non – uniformity** in a specific irradiated area of the sample. The thickness deviation (σ_t/t) is related to the non – uniformity as follows:

$$\frac{\sigma_t}{t} = \frac{\sigma_{non-unif}}{\Delta E} = \frac{\sqrt{\sigma_{meas}^2 - \sigma_{sim}^2}}{\Delta E}$$


M. Fisichella et al., Eur. Phys. J. A 61, 144 (2025)

- Irradiation of different areas of the sample to determine the **global non – uniformity**




Global non – uniformity \equiv deviation thickness between the different irradiation points

Target characterisation at INFN – LNS: characterisation of a multilayer graphene sample




Experimental α – spectrum presented in blue in comparison with the simulated one presented in red. From the σ_{meas} and σ_{sim} we can determine the $\sigma_{\text{non-unif}}$

Typical experimental α – spectra in 5 different irradiated points of the same sample




Target characterisation at INFN – LNS: characterisation of a target evaporated onto a substrate

- 1) APS and characterisation of the substrate
- 2) Evaporation of the target material onto the substrate
- 3) APS using (target + substrate) and characterisation of the target



C (graphene) 0.321 mg/cm^2 (1.59 μm)
 $dE=247\text{keV}$



Te evaporated on a) with a thickness of 0.253 mg/cm^2
 $dE=77\text{keV}$



Targets for high – intense beams

- Experimental campaigns with high – intensity beams are planned at INFN – LNS
- The case of NUMEN experimental campaign:
 - ❖ ^{18}O and ^{20}Ne beams at energies $>15\text{MeV/u}$, beam intensity up to 10^{13} ions/sec → deposition of **several W/cm² in the target**
 - ❖ Isotopically enriched target will be deposited on a **highly thermally conductive foil of graphene**
 - ❖ A **cooling system** for the target will be used
 - ❖ Target materials: Ca, Ge, Se, Zr, Mo, Cd, Pd, Sn, Te, Sm, Pt and many other with a **typical thickness of $\sim 250 - 500 \mu\text{g/cm}^2$**
 - ❖ Both target and substrate should be as uniform as possible to maintain **high energy resolution**
 - ❖ If possible, we should **avoid buffer materials**
- Intense activity on target development and characterisation is ongoing in collaboration with INFN-Torino, INFN-Genova, University of Genova, INFN-LNL and University of Catania

See also talk by Daniela Calvo this afternoon !



*F. Cappuzzello et al.,
Int. J. Mod. Phys. A
36, 2130018 (2021)*



Target production of NUMEN interest: evaporation in standard conditions

Standard evaporation conditions:

1. Low evaporation rate
2. No backing heating
3. No buffer

Standard evaporation conditions

	Tellurium	Germanium	Selenium	Molybdenum
Evaporator	L300 resistive source	L560 resistive source	L300 resistive source	L560 e-beam
Used material	0.4 g	0.6 g	0.5 g	0.6 g
Distance source – backing	210 mm	250 mm	250 mm	200 mm
Heating substrate	NO	NO	NO	300 °C
Buffer	NO	NO	NO	NO
Evaporation rate	0.2 Å/s	0.2 Å/s	1 Å/s	0.3 - 0.4 Å/s

M. Fisichella et al., Eur. Phys. J. A 61, 144 (2025)



Target production of NUMEN interest: some challenging cases

➤ There are several materials on which the application of standard evaporation conditions did not provide satisfactory results.

➤ Depending on our needs, the conditions may be modified accordingly (e.g. the **Mo** case) with satisfactory results


Protocol published in: Eur. Phys. J. A 61, 144 (2025)

➤ One of the special cases is the **Sn** : development of Sn targets at standard or even at non-standard conditions is a difficult task!



Sn (300 μg/cm²) on HOPG @ 150°




@ TRUSTECH, F. Pinna, PhD thesis





@ LNS, M. Fisichella, PhD thesis
(standard evaporation conditions)

Target production : the case of Sn

- One of the special cases is the **Sn** : development of Sn targets at standard or even at non-standard conditions is a difficult task!




- Sn ($250 \mu\text{g}/\text{cm}^2$) + C ($40 \mu\text{g}/\text{cm}^2$)
- Evaporation at standard conditions
- Similar behaviour for all 5 samples of the same evaporation



- Sn ($250 \mu\text{g}/\text{cm}^2$) + Bi (100 \AA) + C ($40 \mu\text{g}/\text{cm}^2$)
- Similar behaviour for all 6 samples of the same evaporation
- Evaporation with and without Bi buffer to be performed with the new UNIVEX 400 evaporator

Summary and future perspectives

- ❑ The development of the appropriate target is essential for the success of nuclear physics experiments
- ❑ The target characterisation:
 - ✓ Provide a tool to examine if the target produced is indeed suitable for the experiment
 - ✓ is of great importance towards the selection of the most suitable manufacturing procedure
- ❑ Target development and characterisation are essential towards the success of large – scale experimental campaigns with high – intensity beams, scheduled in the near future at INFN – LNS
- ❑ A protocol was established for the characterisation of both self – supported and (target + substrate) samples by means of the determination of local and global non – uniformity
 - To be extended in more elements / isotopes
- ❑ New opportunities with the new evaporator recently installed at INFN – LNS



❑ Turbulences in the global market may seriously affect the target production and therefore, the timely schedule of nuclear physics experiments

Targets database within the EURO-LABS project

**Radia Rahali¹, Vasilis Soukeras^{2,3}, Manuela Cavallaro², Christelle Stodel¹, Juan Esposito⁴
for the WP2.5.2 Collaboration**

¹Grand Accélérateur National d'Ions Lourds (GANIL), Caen, France

²Istituto Nazionale di Fisica Nucleare – Laboratori Nazionali del Sud (INFN – LNS), Catania, Italy

³University of Catania, Catania, Italy

⁴Istituto Nazionale di Fisica Nucleare – Laboratori Nazionali di Legnaro (INFN – LNL), Legnaro, Italy

*Second meeting on “Targets for Nuclear Physics”
GSI Helmholtzzentrum für Schwerionenforschung, Darmstadt, Germany*

Targets for Nuclear Physics

□ What are the main properties of a target that we are interested in ?

- ✓ Thickness
- ✓ Homogeneity / Uniformity
- ✓ Purity
- ✓ Isotopic or natural material
- ✓ Possible existence of buffer or substrate material



Strongly depends on the experiment !

Why a new target database is important ?

- 🎯 Find a target not available in – home & facilitate collaboration and target sharing
- 🎯 Improve fabrication techniques
- 🎯 Improve characterisation techniques
- 🎯 Exchange valuable information on target properties, manufacturing and characterisation among different scientific communities that use the same materials.
- 🎯 Collect and centralize data on target production and properties, including both existing and newly developed targets across laboratories.



Final goal: a significant improvement in service to the nuclear physics community - strengthening our knowledge and relationships, and facilitating exchange, collaboration, and knowledge transfer at the European level.



What is the status of the database ?

❖ Database (continuously updated) is already available in excel file including:

- Target material (with Z, A)
- Isotopic or natural material
- Thickness
- Chemical form (e.g. metal)
- Substrate (if exists)
- Manufacturing technique
- Laboratory
- Characterisation technique (if exists)
- Status (already done or in R&D)
- Applications (medical, nuclear physics, industry, material science, ...)
- Publications / references



❖ **Final goal:** database to be publicly accessible to the scientific research community as an online interactive chart



Why It is so important

- Target preparation is a key step toward successful nuclear physics experiments
- A shared database ensures **reproducibility, efficiency, traceability, and quality**
- The laboratories will be able to **contribute, consult, and collaborate** through a digital hub, **expanding** the database

See also talk by
Christelle Stodel
this afternoon !